Description

Light-cure, highly resistant composites for core build-up.

Core D LC is available in a tooth-colored shade which is ideal in situations where esthetics of the core is important, and in a white shade when contrast is needed.

Indications

Core build-up of vital or non-vital tooth.

Safety notes**Warnings**

- Contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

Precautions

- The use of this product is restricted to dental professionals.
- The product must not be used for any purpose other than specified in the indications.
- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindication

Product is contraindicated in patients with a history of allergic reaction to methacrylate resins.

Composition

Hydrophilic dimethacrylates - modified urethane dimethacrylates - barium aluminosilicate glass - amorphous silica - Camphorquinone (CQ) - accelerators - stabilizers - pigments

Incompatibilities

Do not use eugenol-containing materials in the teeth to be restored with Core D LC, since the eugenol may cause discoloration and may interfere with hardening and cause softening of the polymeric components of the materials.

Clinical procedures**Core built up on vital tooth**

1. Core D LC should be used use within 5 minutes from dispensing, to avoid risk of polymerization due to operating light or day light.
2. Prepare the cavity in the usual manners.
3. Clean the tooth structures, rinse thoroughly with a water spray, and lightly air dry. Do not over dry.
4. In case of a cavity close to the pulp, an indirect capping of hard calcium hydroxide should be applied to avoid any post operative sensibility. Do not use eugenol-based cement.
5. Apply an etching gel, following manufacturer instructions if a non self-etching bonding is used.
6. The surfaces must be kept uncontaminated before using a self-etching bonding.
7. The use of a rubber dam is highly recommended to avoid saliva contamination. In case a non-self-etching bonding system is used, the rubber dam should be placed after the surfaces will have been treated.
8. Refer to manufacturer's instructions for tooth conditioning and adhesive application. For optimum bonding, we recommend the use of Healbond MP (self-etching single bottle universal bonding) or Healbond Max MP (self-etching two bottles universal bonding).
9. If needed, use transparent matrix strip and wedges.
10. Remove syringe cap and install an intra-oral tip on the syringe.
11. Gently depress syringe plunger, to begin the flow of material. Do not use excessive force.
12. Apply Core D LC in 2 to 6mm layer maximum, directly into the prepared cavity.
13. Polymerize each layer (from 20 seconds for 2mm, to 40 seconds for 6mm), at a light intensity of at least 800 mW/cm².
14. Just after polymerization is completed, use rotary instruments to finish the crown preparation (final shaping and finishing) in the same clinical trial.
15. After use, leave the intra-oral tip in place on the syringe. It will protect the nozzle until the next application.

Core built up on non-vital tooth

1. Remove temporary fillings from the cavity.
2. Prepare the cavity and shape the root canals in the usual manners. Check that the walls of root canals are free of any materials to insure a proper bonding.
3. Rinse and dry all surfaces.
4. Apply an etching gel, following manufacturer instructions if a non-self-etching bonding is used.
5. The surfaces must be kept uncontaminated before using a self-etching bonding.
6. The use of a rubber dam is highly recommended to avoid saliva contamination. In case a non-self-etching bonding system is used, the rubber dam should be placed after the surfaces will have been treated.

Cementation of endodontic fiber-post

1. The chosen post is fitted into the canal.
2. Treat surface of post according to adhesive manufacturer's and/or dental laboratory's instructions.
3. Cement the post using INNOCEM or any other definitive dual-cure luting cement, according to the IFU.

Core buil up

1. Refer to manufacturer's instructions for tooth conditioning and adhesive application. For optimum bonding, we recommend the use of Healbond MP (self-etching single bottle universal bonding) or Healbond Max MP (self-etching two bottles universal bonding).
2. If needed, use transparent matrix strip and wedges before applying Core D LC.
3. Core D LC should be used use within 5 minutes from dispensing, to avoid risk of polymerization due to operating light or day light.
4. Remove syringe cap and install an intra-oral tip on the syringe.
5. Gently depress syringe plunger, to begin the flow of material. Do not use excessive force.
6. Apply Core D LC in 2 to 6mm layer maximum, directly into the prepared cavity.
7. Polymerize each layer (from 20 seconds for 2mm, to 40 seconds for 6mm), at a light intensity of at least 800 mW/cm².
8. Just after polymerization is completed, use rotary instruments to finish the crown preparation (final shaping and finishing) in the same clinical trial.
9. After use, leave the intra-oral tip in place on the syringe. It will protect the nozzle until the next application.

Storage & shelf-life

- The product must be used by the expiration date indicated on the package.
- Do not expose the product to sun-light or warm sources.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, in order to restore its normal viscosity and curing time.
- Use the product at room temperature.

Important notice

- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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