# Cirus Flow

**Définition**  
Light-cure flowable composite.  
**Safety notes**  
**Warnings**

- Contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

- Precautions**
- The use of this product is restricted to dental professionals.
  - The product must not be used for any purpose other than specified in the indications.
  - Wear suitable protective eyewear, mask, clothing and gloves.
  - Protective eyewear is recommended for patients.

**Contraindication**  
Product is contraindicated in patients with a history of allergic reaction to methacrylate resins.

**Composition**  
Hydrylic dimethacrylates - modified urethane dimethacrylates - barium aluminosilicate glass - amorphous silica - Camphorquinone (CQ) - accelerators - stabilizers - pigments

**Clinical procedures**  
**Restoration of deficient enamel**

1. Prepare the surface of the tooth by prophylaxis, using a non-fluoride prophy-paste. Wash and dry.
  2. Apply etching gel onto the enamel and dentine surfaces beginning with the enamel bevels. Rinse for 20 seconds with water. Dry it in with a water and oil-free air-stream, but do not desiccate. A slightly wet dentine surface is important for the function of last generation bondings. Properly etched tooth should have a dull, chalky-white appearance.
  3. Isolate the tooth to be restored from salivary contamination (rubber dam). Select the desired VITA shade with the tooth slightly moist.
  4. Apply the light-curing bonding agent of your choice. Or use a self-etching bonding agent.
  5. Apply Cirus Flow in layer of 2mm max. over the prepared enamel.
  6. Light cure with a lamp (min. power : 500 mW/cm²) for 20 seconds. Hold the nose of the lamp as close as possible to the surface of the composite.
  7. Use finishing tools as for composites. On heavily discolored teeth ( e.g : tetracycline stains, fluorosis, radiological treatment, etc), the application of an opaque, masking layer is recommended. In such cases, the restorative procedure is modified as follows : After cleaning, etching and bonding, as described above, apply a thin layer of the opaque as an under-layer and cure for 20 seconds. Apply the final coating of Cirus Flow over the opaque layer immediately after cure of the latter. For best adhesion, the opaque layer should be left unfinished or its surface should be roughened.
- Restorations of shallow Class V cavities**
1. Prepare the cavity in a conservative manner, removing minimum tooth structure. Bevel the enamel adjacent to the exposed dentin.
  2. Etch the dentin and enamel, wash, dry and isolate from salivary contamination.
  3. Fill the cavity and finish the restoration as described above.

- Sandwich technique**
1. For MOD and class II, it is recommended to apply a layer of flowable composite prior to fill the cavity with a condensable composite, to get the best possible margin.
  2. Once the tooth preparation, isolation and bonding done, apply a layer of Cirus Flow onto the surface of the cavity, in order to obtain a smooth surface allowing the condensable composite to reach every part of the cavity.
  3. Light cure with a lamp (min. power : 500 mW/cm²) for 20 seconds. Hold the nose of the lamp as close as possible to the surface of the composite.
  4. Then, apply the condensable composite, following manufacturer's instruction for use.

- Pits & fissures obturation**
1. Prepare the surface of the tooth by prophylaxis, using a non-fluoride prophy-paste. Wash and dry.
  2. Apply etching gel onto the enamel. Rinse for 20 seconds with water. Dry it in with a water and oil-free air. Properly etched tooth should have a dull, chalky-white appearance. Highly mineralized tooth may require an additional one to two minutes etching to achieve this effect.
  3. Isolate the tooth to be restored from salivary contamination (rubber dam).
  4. Apply a little quantity of Cirus Flow of the desired VITA shade at the dry surface of the fissure. Using a disposable brush (Fine or XFine size), push the flowable composite down to the bottom of the fissure.
  5. Light cure with a lamp (min. power : 500 mW/cm²) for 20 seconds. Hold the nose of the lamp as close as possible to the surface of the composite.
  6. Finish and polish the restoration material.

**Cirus Flow can also be used in conjunction with a bonding system, to splint mobile teeth with fiber glass.**

- Storage & shelf-life**
- The product must be used by the expiration date indicated on the package.
  - Do not expose the product to sun-light or warm sources.
  - Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, in order to restore its normal viscosity and curing time.
  - Use the product at room temperature.
- Important notice**
- Do not throw away the box and the instruction for use until the product is finished.
  - High humidity could damage the packaging on which important instructions and legal information are printed.
  - Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
  - Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

**Disclaimer**  
These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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