Definition

Dual-cure, self-etch, self-adhesive cement, universal use, for definitive cementation of all types of indirect restorations, whatever their composition. The extra-fine intraoral tips are an answer to all clinical situations.

Innocem can be used either as a "Total Etch" or a "Self-Etch" procedure. The additional etching of the tooth structure with G-Etch increases the adhesive strength.

Indications

- Definitive cementation of ceramic, Zirconia, composite or metal crowns and bridges, on natural teeth or implants.

In combination with Healbond MP, universal bonding, Innocem is indicated for:

- Definitive cementation of metal, ceramic, or composite inlays, onlays and veneers, with a low retention. Etching of dental structure is highly recommended.
- Definitive cementation of Maryland bridges. Etching of dental structure (especially enamel) is necessary.
- Definitive cementation of endodontic metal or fiber posts.

Safety notes**Warnings**

- Contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

Precautions

- The use of this product is restricted to dental professionals.
- The product must not be used for any purpose other than specified in the indications.
- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindication

Product is contraindicated in patients with a history of allergic reaction to methacrylate resins.

Composition

Phosphoric acid monomer MDP - Hydrophilic dimethacrylate - Urethane dimethacrylate (UDMA) - Barium aluminosilicate glass - Camphorquinone (CQ) - Accelerators - Stabilizers - Nanofillers - Pigments

Clinical procedures**First use of a new automix-syringe**

- Remove and discard the cap of the automix-syringe. Under no circumstances should the cap of the automix-syringe be used to reseal the automix-syringe as this could result in the inclusion of air bubbles, inducing a cross-contamination (base/catalyst) and leading to a prematurely polymerization of the product.

- Check whether both automix-syringe openings are free and remove obstructing material, if any, with a suitable instrument.
- Before installing a mixing tip, bleed the automix-syringe onto a mixing pad until the base and catalyst pastes are at the same level.

Tooth structures preparation

1. Clean the tooth structures, rinse with a water spray, and lightly air dry. Do not over dry. The tooth structures should be just dry enough that the surface has a slightly glossy appearance.
2. For veneers and Maryland bridges, tooth structure must be etched.

Treatment of metal, composite, Zirconia and aluminum oxide ceramic restorations

1. Blast or have the composite, metal, Zirconium or non-etchable ceramic restoration surfaces to be luted, blasted with aluminum oxide by the technician.
2. Clean the blasted surfaces with alcohol and dry with air free of water and oil.

Treatment of ceramic restorations

1. If the prosthodontic element has not been prepared by the technician (etching or silanisation), use hydrofluoric acid to etch the glass ceramic restoration surface to be luted. A silanisation "at the chair" could also be done.
2. Rinse thoroughly with water for 15 seconds and dry with air free of water and oil.

Importante notice

1. If the restorations are Zirconia, ceramic or composites onlays, inlays or veneers, with low retention, or Maryland bridges and low abutment teeth, the etching of dental structures is highly recommended (necessary for Maryland bridges).
2. To improve even more the adhesive strength of Innocem, apply Healbond MP, universal multi-purpose adhesive on prepared dental structures and restorations, according to the instructions for use.
3. Wait 15 seconds before drying slightly for 5 to 7 seconds until the solvent has evaporated. Light cure 5 to 10 seconds if you do not want the adhesive "to move" on the surface. This last step is not mandatory.
4. The use of a dual activator is not necessary when Healbond MP is used in conjunction with Innocem, since the cement will turn the adhesive automatically "dual".

Preparation of Root Canals & Treatment of glass fiber-reinforced posts

1. Choose and try the selected post, in the canal.
2. Clean the root canal with a solution of sodium hypochlorite (NaOCl). Rinse immediately and thoroughly with water, dry with paper tips.
3. Clean the post with alcohol and dry it with air free of water and oil. When using Elsdont Fiber Post, application of a silane is not required. Other fiber-reinforced posts should be treated in accordance with their instructions for use.
4. Apply Healbond MP on both the post and canal walls, following IFU and wait 15 seconds before drying slightly for 5/7 seconds until the solvent has evaporated. Do not light-cure.

Cementation

To avoid pulp irritation on vital teeth, cover deep areas, close to the pulp with an anti-bacterial liner [Headlent], or a calcium hydroxide liner.

1. When using a new syringe of Innocem, remove dual syringe cap, dispense and discard a small amount of material. Be sure the material is flowing freely from both ports.
2. Attach a mixing tip on the dual syringe.

Cementation of indirect restorations (inlays, onlays, crowns and bridges)

1. Remove the provisional restoration.
2. Remove temporary cement residus of the prepared tooth or implant abutment surface.

3. Rinse thoroughly.
4. Dry without excess the prepared tooth or implant abutment.
5. Try in restoration to verify its good adaptation, especially at the limits of the preparation.
6. Depress syringe plungers to begin the flow of material. Do not force. Extrude enough product to deposit a thin uniform layer in the inside of the definitive restoration.
7. Seat the restoration carefully onto the prepared tooth or implant abutment, leaving overflow cement.
8. Light-cure the excesses of cement on the edges, for 3 to 5 seconds until it will turn to a gel. Remove them with a probe.
9. For composite and ceramic restorations, light-cure the cement through the restoration.
10. Hold the restoration firmly in place for 2 minutes 30 seconds max. (self-cure mode).
11. Polish the joint with strips or rubber discs cups.
12. Check occlusion.

Cementation of a post in a root canal**Option 1**

1. Use Healbond MP in conjunction with Core D, core build-up material (see IFU of these products).

Option 2

1. Attach an Xfine intra-orally to the mixing tip for direct application in the root canal.
2. Insert the tip as deeply as possible in the root canal and apply Innocem, beginning apically. Keep the intraoral tip in the cement and slowly move upwards.
3. Place the post in the root canal filled with cement and apply moderate pressure to hold the post in position and rotate the post slightly during insertion to avoid the inclusion of air bubbles.
4. Let the cement flowing out of the canal.
5. Light-cure Innocem for 5 seconds to remove the excesses.

Cementation of veneers

1. In cases of low retention, apply Healbond MP on the etched surfaces of the tooth, lightly rubbing the dentin surface for 15 seconds Dry slightly for 5 to 7 seconds until the solvent has evaporated, and light-cure for 5 to 10 seconds if you do not want the adhesive "to move" on the surface. This last step is not mandatory.
2. Inject directly Innocem into the veneer surface. Gently place the veneer on the tooth leaving the cement overflow.
3. Cure through the veneer in place, on the vestibular face, and distance to margins, using the "boost power" of the Led curing light, for about 5 seconds.
4. After having removed cement excess, light-cure the entire surface, as well as the joint, for 10 seconds.

Setting time

- Total setting time is 3 minutes and 10 seconds, (self-cure mode) at 37°C, and 20 seconds (translucent shade) to 30 seconds (opaque shades) in the light-cure mode.

Storage & shelf-life

- The product must be used by the expiration date indicated on the package.
- Do not expose the product to sun-light or warm sources.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, in order to restore its normal viscosity and curing time.
- Use the product at room temperature.

Presentations

Translucent-A2 shade - Opaque A2 shade - Opaque A3 shade

Important notice

- Leave the used mixing tip in place on the automix-syringe until next use.

- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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