# H-Bond 5 ml

MADE IN FRANCE

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GB

**Définition**

H-Bond is a single component self-etching, light-cure universal bonding system.

Depending on the case, this adhesive system can be used:

- In a "self-etching" mode, to enable the shortest possible treatment time and to minimize postoperative sensitivities.
- With selective enamel etching to maximize the adhesion to the enamel.
- With enamel and dentin etching ["total etch"], if the etching process cannot be restricted to the enamel, or when the tooth is devitalized.

**Indications for use**

- To bond light-cure composite during direct restoration procedures.
- Intraoral repair of indirect composite restorations, porcelain fused to metal, and all-ceramic restorations.
- In combination with Universal Activator to modify the curing mechanism from light-curing to dual-curing mode, to bond dual-cure and/or self-cure core build-up materials to the walls of cavities and to the surface of metal posts and reinforced fiber posts.
- In combination with Universal Activator and Innocem (or any other resin-based definitive cement) for:
  - Definitive cementation of small Maryland bridges (2 to 3 elements).
  - Definitive cementation of composite and all-ceramic veneers when the retention is poor.
  - Definitive cementation of metal posts or reinforced fiber posts.
- Protective varnish for glass ionomer fillings.

**Safety notes**

**Warnings**

- The product contains methacrylate resins and alcohol which may irritate skin, eyes and oral mucosa.
- It may cause damage to the pulp and allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

- The product is a corrosive liquid as it contains ethanol and is flammable. Keep away tram sources of ignition.

**Précautions**

- The use of this product is restricted to dental professionals.
- The product must not be used for any purpose other than specified in the indications.
- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

**Contraindication**

H-BOND is contraindicated in patients with a history of allergic reaction to methacrylate resins.

**Device main components**

Silane - phosphoric acid monomer (MDP) - hydrophilic dimethacrylate - Dimethacrylate resins - camphorquinone (CQ) - accelerators - stabilizers - nanofillers - ethanol and water.

**Clinical procedures**

**DIRECT RESTORATIONS**

**Bonding Composites to dentin/enamel**

- Rubber dam is the recommended method if isolation to keep the working area free of blood and saliva, and to protect soft tissues.
- Prepare the cavity with minimal tooth reduction. Cavity floor very close to pulp should be covered with a thin layer
- of hard setting calcium hydroxide material or with an anti-microbial liner (e.g: Healdent).
- Rinse thoroughly with water spray and gently air dry (do not over-dry).
- The tooth structures should be just dry enough that the surface has a slightly glossy appearance.
- Apply a uniform coating of H-Bond onto the prepared dentine and enamel surfaces for 15 seconds with a rubbing motion.
- Gently dry with water-free and oil-free air for 10 to 20 seconds to evaporate the solvent base and to disperse the bonding homogeneously. Keep the surface slightly wet for maximum bonding strength.
- Light cure the bonding for 10 seconds.
- Apply composite of your choice per manufacturer’s instructions for use.

**Note**

H-BOND produces high enamel bond strengths without the use of etching. However, etching gel can be used on enamel and/or dentin if a total etch or selective enamel etch technique is desired. In that case, during step 2, after preparation and before applying a cavity liner, if needed, etch enamel or dentin surfaces for 15 seconds, then rinse thoroughly with water spray for 15 seconds and gently air dry (do not over-dry). Proceed to step 3.

**INDIRECT RESTORATIONS with poor retention (Veneers, Maryland bridge)**

**Important notice**

- In case of indirect restorations with normal retention, the use of HBond, as a complement of Innocem, is not necessary.

- Light-curing of the mixture (HBond + Universal Activator) is optional if Innocem resin-cement is used for cementation.

**Tooth structures preparation**

- Prepare the teeth in a conventional manner. Enamel surfaces should be etched, to improve the adhesive strength.
- Try in restoration to verify its good adaptation, especially at the limits of the preparation.
- Clean the tooth structures, rinse thoroughly with a water spray, and gently air dry. Do not over dry.
- The tooth structures should be just dry enough that the surface has a slightly glossy appearance.
- Mix a drop of H-Bond adhesive and a drop of Universal Activator for 5 seconds.
- Apply a uniform coating of the mixture onto the prepared dentin and enamel, using a disposable applicator, for 15 seconds with a rubbing motion.
- Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
- Light cure for 5 to 10 seconds. That will also stop the mixture "moving" on the surface.

**Preparation of restorations**

**Treatment of Lithium Disilicate, Ceramic/Porcelain veneers [etchable]**

- Abrase the internal surface of the restorations with 50µm alumina with a pressure of about 30 psi (0.2 MPa),
- Etch the restoration with hydrofluoric acid for 1 minute.
- Rinse thoroughly with water for 15 seconds and dry with air, free of water and oil.
- Mix a drop of H-Bond adhesive and a drop of Universal Activator for 5 seconds.
- Apply a uniform coating of the mixture on the internal surface of the prepared restorations for 15 seconds with a rubbing motion.
- Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
- Light cure for 5 to 10 seconds. That will also stop the mixture "moving" on the surface.

**Treatment of composite veneers [non-etchable]**

- Abrase the internal surface on the restorations with 50µm alumina with a pressure of about 15 psi (0.1 MPa).
- Clean the prepared surfaces with alcohol and dry with a water-free and oil-free air.
- Mix a drop of H-Bond adhesive and a drop of Universal Activator for 5 seconds.
- Apply a uniform coating of the mixture on the internal surface of the prepared restorations for 15 seconds with a rubbing motion.
- Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
- Light cure for 5 to 10 seconds. That will also stop the mixture "moving" on the surface.

**Treatment of Maryland Bridges**

- Abrase the internal surface of the bridge with 50µm alumina with a pressure of about 60 psi (0.4 MPa).
- Clean the prepared surfaces with alcohol and dry with a water-free and oil-free air.
- Mix a drop of H-Bond adhesive and a drop of Universal Activator for 5 seconds.
- Apply a uniform coating of the mixture on the internal surface of the prepared restorations for 15 seconds with a rubbing motion.
- Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
- Light cure for 5 to 10 seconds. That will also stop the mixture "moving" on the surface.

**Cementation and seating of restoration**

**Veneers, Maryland bridges**

- Rubber dam is the recommended method if isolation to keep the working area free of blood and saliva, and to protect soft tissues.
- Apply Innocem, or any other resin-based definitive cement of your choice, to the restoration and seat it.
- Cure the resin-based definitive cement per manufacturer’s instructions for use.

**Note**

Ceramic restorations must be cemented before checking occlusion (risk of breakage)

**Definitive cementation of root canal posts and bonding of dual-cure core build up composite**

**Preparation of root canals and posts**

- Rubber dam is the recommended method if isolation to keep the working area free of blood and saliva, and to protect soft tissues.
- Clean the prepared root canal with a 5% solution of sodium hypochlorite (NaOCl).
- Rinse thoroughly with water spray. Dry with paper tips.
- Trial fit the post. It should fill 2/3 of the root canal.
- Clean the post with alcohol and dry it with water-free and oil-free air.
- Mix a drop of H-Bond and a drop of Universal Activator, for 5 seconds.
- Apply the mixture on the walls of the root canals for 15 seconds with a rubbing motion.
- Apply a thin layer of the mixture on the post.
- Slightly dry both the walls and the post for 10 seconds to evaporate the solvent base and to disperse the bonding homogeneously. Keep the surface slightly wet for maximum bonding strength.
- Light cure both the surface of the walls and the post for 5 to 10 seconds.
- Using a dual cure core built up composite of your choice, cement the post and fill the whole cavity, per manufacturer’s instructions.

**Note**

Innocem or any other resin-based definitive cement can be used, per manufacturer’s instructions, to cement the

post. The core built up material will be used to fill the cavity only.

**Intra-oral repair of indirect composite, fused-ceramic to metal and all-ceramic restorations**

**Indirect composite restoration**

- Roughen the surface to be repaired (sandblasting) and clean with alcohol. Dry it with water-free and oil-free air.
- Apply a uniform coating of H-Bond onto the prepared surfaces for 15 seconds with a rubbing motion.
- Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the bonding homogeneously. Keep the surface slightly wet for maximum bonding strength.
- Light cure for 10 seconds.
- Apply the composite of your choice, per manufacturer’s instructions.

**Fused-ceramic to metal and all-ceramic restoration**

- Sandblast the metal if necessary and clean with alcohol. Dry it with water-free and oil-free air.
- Etch the restoration with hydrofluoric acid for 1 minute.
- Rinse thoroughly with water and dry with water-free and oil-free air.
- Apply a uniform coating of H-Bond onto the prepared surfaces for 15 seconds with a rubbing motion.
- Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the bonding homogeneously. Keep the surface slightly wet for maximum bonding strength.
- Light cure for 10 seconds.
- For fused ceramic to metal restoration, when metal is visible, apply an opaquer onto the bonding and light-cure for 20 seconds, per manufacturer’s instructions.
- Apply the composite of your choice, per manufacturer’s instructions.

**Storage and shelf-life**

- The product must be used by the expiration date indicated on the package.
- The product should be tightly closed immediately after use.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, in order to restore its normal viscosity and curing time.
- Do not expose the product to sun-light or warm sources.
- Never store close to sparkles and open fire.

**Important notice**

- Do not throw away the box and the instruction for use until the product is completely finished.
- High humidity could damage the packaging on which important instructions and legal information is printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

**Disclaimer**

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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