## FP Root canal fiber posts

GB

### Description

FP Fiber Posts are cylindrical posts with a conical end for the last 2 mm, to give strong resistance and retention. FP Fiber Posts contain an organic matrix (40%) and glass fibers (60%).

### Safety notes

#### Warnings

- FP Fiber Posts contain an organic matrix (40%) and glass fibers (60%) which may irritate eyes and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water, remove contact lenses, and seek medical attention.

#### Precautions

- Use of rubber dam or adequate isolation is recommended.
- Wear suitable protective eyewear, mask, and gloves.
- Protective eyewear is recommended for patients.

### Contraindication

FP Fiber Posts is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components.

### Bruxism

#### Deep overbite

#### Instructions for use

Root canal preparation, choice and luting of FP posts

- Use a X-ray to determine the appropriate diameter and depth of the planned post space preparation. The diameter of the post should be smaller than the diameter of the prepared root canal. For perfect retention, it is important that 2/3 of the length of the post is in the root canal.
- Prepare the root canal by using the drill corresponding to the chosen post (a colour code facilitates the choice of the correct drill).
- Rinse the prepared root canal in order to eliminate debris. Dry with paper points. Try the chosen post in the root canal to check that nothing blocks the post, preventing it from reaching the limit of the preparation.
- Isolate the operating field by using a rubber dam.
- Desinfect the prepared root canal (5% NaOCl) and dry with paper points.
- Depending on the type of bonding used, etch or not the root canal, the dentin wall and the remaining enamel, following manufacturer's instructions. For optimum bonding, we recommend use of HEAL-BOND DUO (with etching) or HEALBOND DUO SE (self-etching).
- Apply the bonding of your choice on a slightly wet surface. Dry the excess of the bonding agent from the canal with paper points.
- Do not polymerize the bonding.
- If needed, adapt a matrix band around the tooth (not necessary when using CORE D SC, CORE D or CORE D FLOW).
- Silanise the chosen post following manufacturer's instructions.
- Apply some CORE D SC (self-cure core build-up material), CORE D (dual-cure core build-up material) or CORE D FLOW (Hi Tech dual-curing core built up material, re-inforced with Zirconium), on the conical end of the post. Do not inject in the canal, to prevent the core built up material mixing with the bonding agent.
- Put the post in place in the canal by pushing it down gently in the desired position. Hold it with a slight pressure.
- For the self-curing composites (CORE D SC), the setting time is 1 minute maximum, beginning of injection.

- For dual-cure composites (CORE D, CORE D FLOW), light-cure for 60 seconds using a curing unit with a light intensity of more than 400 mW/cm<sup>2</sup>. The light tip of the lamp should be positioned as close to the post as possible.

#### Storage & shelf-life

Store in a dry, cool and well-ventilated place, away from humidity.

#### References & presentation

KFP-15 - Intro kit of 3 x 5 posts (1.00 mm – 1.20 mm - 1.35 mm), under individual blister + 3 assorted drills.

#### Important notice

- Do not throw away the box and the instruction for use until the product is completely finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.

#### Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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### Bruxisme

#### Superocclusion incisive/surocclusion

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