HEALBOND MP

GB

Description of the product

HEALBOND MP is a single component self-etching, light-cure universal bonding system.

Depending on the case, it can be used:

- In a "self-etching" mode, to enable the shortest possible treatment time and to minimize postoperative sensitivities.
- With selective enamel etching to maximize the adhesion to the enamel.
- With enamel and dentin etching ("total etch"), if the etching process cannot be restricted to the enamel, or when the tooth is devitalized.

Indications for use

1. To bond light-cure composite during direct restoration procedures.
2. Intraoral repair of indirect composite restorations, porcelain fused to metal, and all-ceramic restorations without primer.
3. To bond Elsodont dual-cure and/or self-cure core build-up materials without using Universal Activator to the walls of cavities and to the surface of metal posts and reinforced fiber posts.
4. In combination with InnoCem, definitive resin cement, as an adhesive for:
 - Definitive cementation of small Maryland bridges [2 to 3 elements].
 - Definitive cementation of composite and all-ceramic veneers when the retention is poor.
 - Definitive cementation of metal posts or reinforced fiber posts.
5. Protective varnish for glass ionomer fillings.

Safety notes

Warnings

- The product contains methacrylate resins and alcohol which may irritate skin, eyes and oral mucosa. It may cause damage to the pulp and allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.
- The product is a corrosive liquid and is flammable. Keep away from sources of ignition.
- Avoid ingestion of the product. If signs of digestive disorder appear, contact a doctor.

Precautions

- The use of this product is solely restricted to dental personnel.
- The product must not be used for any purpose other than specified in the indications.
- Dental professional should wear suitable protective eyewear, mask, clothing, and gloves.
- Protective eyewear is recommended for patients.
- Isolation of the operating field is requested [e.g. dental dam].
- Dispose according to the national and/or federal standards in force.

Contraindication

The product is contraindicated in patients with a history of allergic reaction to any of the components.

Importante notice

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Incompatibilities

- The product should not be used with eugenol-containing materials for pulp protection or temporary sealing since the eugenol can retard or impair the curing process.
- The product cannot be used with a hydrogen peroxide solution for cleaning cavities since it may weaken the bond strength to the tooth structure.
- The product cannot be used with hemostatic agents containing ferric compounds, since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva, due to remaining ferric ions.

Composition

Silane [0,5-2%] - MDP [5-20%] - Dimethacrylate resins [15-25%] - Camphorquinone [<1%] - Accelerators [>1%] - BHT [<1%] - Glass fillers [1-10%] - Ethanol [15-30%] - Water [15-30%]

Times

Halogen/LED lamp [800 to 1400mW/cm²]: 10 seconds for polymerization. 5 seconds when use with compatible products such as InnoCem or Core D which complete the curing process of HEALBOND MP.

Clinical procedures

DIRECT RESTORATIONS

Bonding Composites to dentin/enamel

1. Rubber dam is the recommended method if isolation to keep the working area free of blood and saliva, and to protect soft tissues.
2. Prepare the cavity with minimal tooth reduction. Cavity floor very close to pulp should be covered with the protective liner of your choice.
3. Rinse thoroughly with water spray and gently air dry (do not over-dry). The tooth structures should be just dry enough that the surface has slightly glossy appearance.
4. Apply a uniform coating of HEALBOND MP onto the prepared dentin and enamel surfaces for 15 seconds with a rubbing motion.
5. Gently dry with water-free and oil-free air for 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
6. Light cure for 10 seconds.
7. Apply composite of your choice per manufacturer's instructions for use.

Note

HEALBOND MP produces high enamel bond strengths without the use of etching. However, etching gel can be used on enamel and/or dentin if a total etch or selective enamel etch technique is desired. In that case, during step 2, after preparation and before applying a cavity liner, if needed, etch Enamel or Dentin surfaces for 15 seconds, then rinse thoroughly with water spray for 5 seconds and gently air dry (do not over-dry).

INDIRECT RESTORATIONS

Note

Light-curing of HEALBOND MP is optional if INNOCEM resin-cement is used for cementation.

Tooth structures preparation

1. Prepare the teeth in a conventional manner. Enamel surfaces should be etched, to improve the adhesive strength.
2. Try in restoration to verify its good adaptation, especially at the limits of the preparation.
3. Clean the tooth structures, rinse thoroughly with a water spray, and gently air dry. Do not over dry. The tooth structures should be just dry enough that the surface has a slightly glossy appearance.
4. Isolate the operating field (dental dam).
5. Apply a uniform coating of HEALBOND MP onto the prepared dentin and enamel, using a disposable applicator for 15 seconds with a rubbing motion.

6. Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
7. Slightly dry both the walls and the post for 10 seconds to evaporate the solvent base and to disperse the bonding homogeneously. Keep the surface slightly wet for maximum bonding strength.

Note

HEALBOND MP produces high enamel bond strengths without the use of etching. However, etching gel can be used on enamel and/or dentin, especially in case of veneers with poor retention, and Maryland bridges. In that case, after step 1.3, etch Enamel or Dentin surfaces for 15 seconds, then rinse thoroughly ensuring that all etchant has been removed. Proceed to step 4.

Preparation of restorations

Treatment of Lithium Disilicate, Ceramic/Porcelain veneers [etchable]

1. Abrase the internal surface of the restorations with 50μm alumina with a pressure of about 30 psi [0.2 MPa].
2. Etch the restoration with hydrofluoric acid for 1 minute.
3. Rinses thoroughly with water for 15 seconds and dry with air, free of water and oil.
4. Apply a uniform coating of HEALBOND MP on the internal surface of the prepared restorations for 15 seconds with a rubbing motion.
5. Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
6. Light cure for 10 seconds.
7. Apply the composite of your choice, per manufacturer's instructions.

Note

HEALBOND MP can be used as silane replacement. Treatment of Metal-based [Maryland Bridge], Zirconia-based, Alumina-based and Composite Restorations [non-etchable]

1. Abrase the internal surface of the restorations with 50μm alumina with a pressure of about 15 psi [0.1 MPa].
2. Clean the prepared surfaces with alcohol and dry with a water-free and oil-free air.
3. Apply a uniform coat of HEALBOND MP on the internal surface of the prepared restoration for 15 seconds with a rubbing motion.
4. Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
5. Light-cure for 10 seconds.
6. For fused ceramic to metal restoration, when metal is visible, apply an opaquer onto the bonding and light-cure it for 20 seconds, per manufacturer's instructions.
7. Apply the composite of your choice, per manufacturer's instructions.

Cementation and seating of restoration

Veneers, Crowns, bridges, inlays, onlays and metal-based restorations [Maryland bridges]

1. Rubber dam is the recommended method if isolation to keep the working area free of blood and saliva, and to protect soft tissues.
2. Apply InnoCem, or any other resin cement of your choice to the restoration and seat it.
3. Cure the cement per manufacturer's instructions for use.

Note

Ceramic restorations must be cemented before checking occlusion [risk of breakage].

Definitive cementation of root canal posts and bonding of dual-cure core build up composite

Preparation of root canals and posts

1. Rubber dam is the recommended method if isolation to keep the working area free of blood and saliva, and to protect soft tissues.
2. Clean the prepared root canal with a 5% solution of sodium hypochlorite (NaOCl).
3. Rinse thoroughly with water spray. Dry with paper tips.

4. Trial fit the post. It should fill 2/3 of the root canal.
5. Clean the post with alcohol and dry it with water-free and oil-free air.
6. Apply HEALBOND MP on the walls of the root canals for 15 seconds with a rubbing motion.
7. Also apply a thin layer of HEALBOND MP on the post.

8. Slightly dry both the walls and the post for 10 seconds to evaporate the solvent base and to disperse the bonding homogeneously. Keep the surface slightly wet for maximum bonding strength.
9. Apply InnoCem, or any other resin cement of your choice. Apply CORE D orany dual cure core-build up material of your choice onto the post, seat the post, vibrate the post slightly avoiding the possibility of trapped air, and then cure the resin cement or the dual cure composite per the manufacturer's instructions
10. Proceed with a core build-up material of your choice in accordance with manufacturer's instructions.

11. Abrase the internal surface of the restorations with 50μm alumina with a pressure of about 30 psi [0.2 MPa].
12. Etch the restoration with hydrofluoric acid for 1 minute.
13. Rinses thoroughly with water for 15 seconds and dry with air, free of water and oil.
14. Apply a uniform coating of HEALBOND MP on the internal surface of the prepared restorations for 15 seconds with a rubbing motion.
15. Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
16. Light cure for 10 seconds.
17. Apply the composite of your choice, per manufacturer's instructions.

Indirect composite restoration

Treatment of indirect composite, fused-ceramic to metal and all-ceramic restorations

Intra-oral repair of indirect composite, fused-ceramic to metal and all-ceramic restorations

Indirect composite restoration

1. Roughen the surface to be repaired [sandblasting] and clean with alcohol. Dry it with water-free and oil-free air.
2. Apply a uniform coating of HEALBOND MP onto the prepared surfaces for 15 seconds with a rubbing motion.
3. Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Keep the surface slightly wet for maximum bonding strength.
4. Light cure for 10 seconds.
5. Apply the composite of your choice, per manufacturer's instructions.

Fused-ceramic to metal and all-ceramic restoration

Indirect composite restoration

1. Sandblast the metal if necessary and clean with alcohol. Dry it with water-free and oil-free air.
2. Etch the restoration with hydrofluoric acid for 1 minute.
3. Rinse thoroughly with water and dry with water-free and oil-free air.
4. Apply a uniform coating of HEALBOND MP onto the prepared surfaces for 15 seconds with a rubbing motion.
5. Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously.
6. Light-cure for 10 seconds.
7. Apply the composite of your choice, per manufacturer's instructions.

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Storage and shelf-life

- The product must be used by the expiration date indicated on the package.
- Do not expose the product to sun-light or warm sources.
- The product should be tightly closed immediately after use.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, to restore its normal viscosity and curing time.
- Use at room temperature.
- Never store close to sparkles and open fire.

Important notice

- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information is printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

Disposal must be carried out in accordance with local, regional and national regulations in force. These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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