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GB

Definition

Light-cure resin-based liner – Radiopaque - HI-SEAL chemically bonds to adhesive primers, composites and other resin-based materials and it micromechanically bonds to dentin.

Safety notes

Warnings

- Contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

Precautions

- The use of this product is restricted to dental professionals.
- The product must not be used for any purpose other than specified in the indications.
- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.

Contraindication

Product is contraindicated in patients with a history of allergic reaction to methacrylate resins.

Composition

Hydrophilic dimethacrylates - modified urethane dimethacrylates - Fluoroaluminosilicate - barium aluminosilicate glass - amorphous silica - Camphorquinone (CQ) - accelerators - stabilizers - pigments

Clinical procedures

1. Prepare, rinse and dry the cavity.
2. Isolate treatment area from salivary contamination (cotton rolls, rubber dam etc...)
3. Apply the material directly over the dentin through a disposable applicator. When increased bond strength to dentin is desirable, apply first your preferred dental adhesive to the dentin. The first coating should not be more than 0.3 - 0.5 mm thick. Remove any material inadvertently applied over enamel margins. Cure for 20 - 25 seconds, keeping the tip as close to the cured surfaces as possible.
4. Apply a second layer, if indicated. The second layer can be up to 1.5mm thick. Cure for 30 seconds.
5. Restore the tooth with a composite, per manufacturer instructions. HI-SEAL is compatible with all commercially available dental restorative materials.

Note

If using HI-SEAL under water-based cements such as glass ionomer, remove any uncured resin from the Hi-Seal surface by wiping it with a pledge of cotton moistened with alcohol.

Storage & shelf-life

- The product must be used by the expiration date indicated on the package.
- Do not expose the product to sun-light or warm sources.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, in order to restore its normal viscosity and curing time.
- Use the product at room temperature.

Important notice

- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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