STILL BITE

Description of the product

Extremely hard, addition-curing bite registration material based on polyvinyl siloxane.

Safety notes

- Warnings**
- Contains polymerizable monomers and surfactants which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
 - Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs immediately remove material with cotton and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.
 - Avoid ingestion/swallowing of material to prevent irritation or obstruction. If accidental swallowing occurs, drink lots of water. This material is not hazardous when small quantities are ingested. Larger quantities may cause bowel obstruction. Seek medical attention in the event of digestive irregularities.

Precautions

- The use of this product is restricted to dental professionals.
- The product must not be used for any purpose other than specified in the indications.
- Wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.
- Higher temperatures reduce working times and setting times (faster), lower temperatures increase them (slower).
- To obtain optimum physical properties, the mix should be completely homogeneous.
- Should be used with properly fitted, rigid trays for most applications. Trays with retentive features (perforations), are recommended.
- Material should extrude easily: do not use excessive force. Excessive pressure may result in unanticipated extrusion of the material or cause cartridge rupture.
- Dispose according to the national and/or federal standards in force.

Contraindication

Product is contraindicated in patients with a history of allergic reaction to any of the components.

Important notice

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Composition

Polydimethylsiloxane with vinyl group, Cristobalite, Pigment

Clinical procedures

First use of a new cartridge

- Remove and discard the cap of the cartridge. Under no circumstances should the cap of the cartridge be used to reseal the cartridge as this could result in the inclusion of air bubbles, inducing a cross-contamination (base/catalyst) and leading to a prematured polymerization of the product.
- Check whether both cartridge openings are free and remove obstructing material, if any, with a suitable instrument.
- Before installing a mixing tip, bleed the cartridge onto a mixing pad until the base and catalyst

pastes are at the same level

Bite registration

1. Remove and discard the used mixing tips.
2. Attach a new mixing tip (yellow ring type) and an intraoral tip (yellow), on the cartridge.
3. Place the cartridge into the dispensing gun.
4. Slightly dry teeth before registration but, if necessary, isolate the surface to be registered with a vaseline film or equivalent (in presence of new composite, for instance).
5. Block important undercuts (notches, embrasures) with wax or liquid dam (ex: Notre Dam).
6. Apply the bite registration material onto dry teeth. The bite registration should be about 3 to 4 cm length.
7. Ask the patient to close light, to obtain a perfect registra-tion of the occlusion.
8. As soon as the material started its polymerization (after about 30 seconds, it reaches an elastic phase), ask the patient to open his mouth.
9. Take the bite registration material off the mouth at this stage. Complete polymerization occurs after 2 more minutes.
10. Trim the bite to keep the cusp tips, with a resin bur, to avoid repositioning mistakes.
11. The bite registration can be disinfected, and its inhibition layer can be removed with alcohol.

Note

This device is used on young adults to elderly people.

Storage & shelf-life

- The product must be used by the expiration date indicated on the package.
- Do not expose the product to sun-light or warm sources.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, in order to restore its normal viscosity and curing time.
- Use the product at room temperature.

Important notice

- Leave the used mixing tip in place on the cartridge until next use.
- Do not touch the material with latex gloves. They can inhibit its polymerization.
- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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