CORE D LC

GB

Description

Light-cure composite for core build-up of devitalized teeth.

Safety notes

Warnings

- Contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.
- Avoid ingestion of the product. If signs of digestive disorder appear, contact a doctor.

Precautions

- The use of this product is restricted to dental personnel.
- The product must not be used for any purpose other than specified in the indications.
- Wear suitable protective eyewear, mask, clothing and gloves for dental professionals.
- Protective eyewear is recommended for patients.
- Isolation of the operating field is requested (e.g. dental dam).
- Dispose according to the national and/or federal standards in force.

Contraindication

Product is contraindicated in patients with a history of allergic reaction to any of the components.

Important notice

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Incompatibilities

The product should not be used with eugenol-containing materials for pulp protection or temporary sealing since the eugenol can retard the curing process.

Composition

Hydrophilic dimethacrylates - modified urethane dimethacrylates (20-30%) - barium aluminosilicate glass (50-60%) - amorphous silica - Camphorquinone (CQ) (<1%) - accelerators - stabilizers - pigments

Clinical procedures

1. Remove temporary fillings from the cavity.
2. Prepare the cavity and shape the root canals in the usual manners. Check that the walls of root canals are free of any materials to insure a proper bonding.
3. Rinse and dry all surfaces.
4. Apply an etching gel, following manufacturer instructions if a non-self-etching bonding is used.
5. The surfaces must be kept uncontaminated before using a self-etching bonding. The use of a rubber dam is highly recommended to avoid saliva contamination. In case a non-self-etching bonding system is used, the rubber dam should be placed after the surfaces will have been treated.

Cementation of endodontic fiber-post

1. The chosen post is fitted into the canal.
2. Treat surface of post according to adhesive manufacturer's and/or dental laboratory's instructions.
3. Cement the post using INNOCEM or any other definitive dual-cure luting cement, according to the IFU.

Core built up

1. Refer to manufacturer's instructions for tooth conditioning and adhesive application. For optimum bonding, we recommend the use of Healbond MP (self-etching single bottle 8th generation, universal bonding) or Healbond Max MP (self-etching two bottles universal bonding).
2. If needed, use transparent matrix strip and wedges before applying Core D LC.
3. Core D LC should be used use within 5 minutes from dispensing, to avoid risk of polymerization due to operating light or day light.
4. Remove syringe cap and install an intra-oral tip on the syringe.
5. Gently depress syringe plunger, to begin the flow of material. Do not use excessive force.
6. Apply Core D LC in 4mm layer maximum, directly into the prepared cavity.
7. Polymerize each layer 40 seconds, using a LED light (> 1000 mW/cm²).
8. Once polymerized, use rotary instruments to complete the crown preparation (final shaping and finishing) in the same clinical session.
9. Put in place a temporary crown until next appointment.
10. After use, leave the intra-oral tip in place on the syringe. It will protect the nozzle until the next application.

Note

This device is used on young adults to elderly people

Storage & shelf-life

- The product must be used by the expiration date indicated on the package.
- Do not expose the product to sun-light or warm sources.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, to restore its normal viscosity and curing time.
- Use the product at room temperature.

Important notice

- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

Dispose of the product in accordance with applicable local, regional and national regulations. This material is intended exclusively for use in dentistry. Processing must be done strictly in accordance with the instructions for use. No liability can be accepted for damage caused by failure to comply with the instructions for use or the prescribed area of application. It is the responsibility of the user to check the material for its suitability and use for a purpose not expressly stated in the instructions for use. Descriptions and information do not constitute a guarantee of properties and are not binding.

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