Phoenix MD

EN

Description of the product

Self-cure resin for temporary crowns and bridges material, with high mechanical properties. Thixotropic, to also make temporary crowns directly in mouth, on the abutment teeth, without a pre-operative impression.

Safety notes

Warnings

- Contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
- Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.
- Avoid ingestion of the product. If signs of digestive disorder appear, contact a doctor.

Precautions

- The use of this product is restricted to dental personnel.
- The product must not be used for any purpose other than specified in the indications.
- Dental professional should wear suitable protective eyewear, mask, clothing and gloves.
- Protective eyewear is recommended for patients.
- Dispose according to the national and/or federal standards in force.

Contraindication

This product is contraindicated in patients with a history of allergic reaction to any of the components.

Important notice

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Incompatibilities

The use on eugenol-containing products with based resin-based dental materials should be avoided since they interfere with the curing process.

Composition

Glass fillers [30-40%] - Dimethacrylate resins [40-60%] - Catalysts <1% - BHT <0.1%] - Pigments <1%]

Times

Total time (working time + setting time) at 22°C, is 2'35" to 2'45".

Clinical procedures

First use of a new cartridge

- Remove and discard the cap of the cartridge. Under no circumstances should the cap of the cartridge be used to reseal the cartridge as this could result in the inclusion of air bubbles, inducing a cross-contamination (base/ catalyst) and leading to a premature polymerization of the product.
- Check whether both cartridge openings are free and remove obstructing material, if any, with a suitable instrument.
- Before installing a mixing tip, bleed the cartridge onto a mixing pad until the base and catalyst pastes are at the same level.

Temporary C&B - Undirect method

- Take a pre-operative impression. Interdental spaces should be rid of silicone material.
- Remove and discard the used mixing tips and attach a new one on the cartridge. Insert the cartridge into the dispensing gun.
- Blow dry impression and fill the concerned teeth area, with Phoenix MD. Always fill from the bottom up, until the impression is ¾ full.
- It is recommended to also inject some resin into the sulcus of the concerned teeth, with a yellow intra-oral tip, to avoid any lack of product, at the limits of the temporary crown.

- Place the impression back in patient's mouth and allow Phoenix MD to set undisturbed for about 40/45 seconds.
- Remove the impression tray from mouth. At that stage, the material will be somehow flexible and will pull over any undercuts or non-parallel preparations. The material will fully cure in 2 additional minutes.
- Remove the temporary crown or bridge from the impression and trim excesses.
- Remove the smear layer from the surface by whipping with a cloth moistened with alcohol.
- Polish and cement the temporary crown with a temporary cement (e.g.: Seal Temp or Seal Temp S).

Temporary C&B - Direct method

- Compared with other C&B temporary resins in cartridges, PHOENIX MD does not flow. The material can be directly injected onto the abutment tooth, to make a perfectly adjusted temporary crown in about 3 minutes maximum.
- Inject PHOENIX MD directly into the sulcus of the abutment tooth that you keep it slightly wet with saliva.
 - Once the tooth is completely covered, ask the patient to bit.
 - After 30 seconds, the material reaches an elastic phase. Take the temporary crown of the abutment tooth and put it back few times, to avoid the temporary crown to be blocked in the inter-proximal undercuts.
 - Shape the temporary crown in mouth with diamond bur, under water spray.
 - Trim the excesses and cement the crown with a temporary luting cement (e.g.: SEALTEMP or SEALTEMP S).

Note

This device is used on young adults to elderly people, since definitive prosthodontic sessions are only conducted on adults.

Storage & shelf-life

- The product must be used by the expiration date indicated on the package.
- Do not expose the product to sun-light or warm sources.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, to restore its normal viscosity and curing time.
- Use the product at room temperature.

Shades

PHOENIX MD is available in the following shades: A1 - A2 - A3 Shades A1, A2 and A3 are based on Vita® shades* *Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany

Important notice

- Leave the used mixing tip in place on the cartridge until next use.
- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

Disposal of the product must be carried out in accordance with local, national and regional regulations in force. These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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