INNOCEM Veneer LC

GB

Description of the product
INNOCEM Veneer LC is a light-cure, radio-opaque heavy flowable nano-hybrid composite-based cement paste. BIO+ formulation BPA-free (no BisGMA/BisEMA), TEGDMA-free and no addition of HEMA. INNOCEM Veneer LC exhibits superior color stability, and it is available in 3 shades: Light [A1], Medium [A2], Bright [Bleach].

The INNOCEM Veneer LC Paste syringe is generally delivered with a Try In Paste syringe of the same shade. INNOCEM Veneer-Try-In Paste is a trial paste to confirm suitable paste shade before placing the restoration. It is water-soluble and can be easily removed after confirming the shade.

Classification: Type2/Class2/Group1, according to EN ISO 4049

Indications
Permanent adhesive luting cement for inlays, onlays and veneers, made of ceramic or composite, with a thickness < 2.0mm and a translucency that enables the complete light-curing of the cement.

Safety notes

- Warnings**
- The product contains methacrylate resins which may irritate skin, eyes and oral mucosa. It may cause allergic contact dermatitis in susceptible persons.
 - Avoid eye contact and prevention irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin.
 - If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.
 - Avoid ingestion of the product. If signs of digestive disorder appear, contact a doctor.

- Precautions**
- The use of this product is solely restricted to dental personnel.
 - The product must not be used for any purpose other than specified in the indications.
 - Dental professional should wear suitable protective eyewear, mask, clothing, and gloves.
 - Protective eyewear is recommended for patients.
 - Isolation of the operating field is requested (e.g. dental dam).

Contraindication
Product is contraindicated in patients with a history of allergic reaction to any of the components.

Important notice
Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

- Incompatibilities**
- The product should not be used with eugenol-containing materials for pulp protection or temporary sealing since the eugenol can retard or impair the curing process.
 - The product cannot be used with a hydrogen peroxide solution for cleaning cavities since it may weaken the bond strength to the tooth structure.
 - The product cannot be used with hemostatic agents containing ferric compounds, since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva, due to remaining ferric ions.
 - Do not use a sodium hypochlorite solution since it may weaken the bonding to the tooth structures.

Composition
Hydrophobic dimethacrylate [10-20%] - Hydrophitic dimethacrylates [1-5%] - Modified urethane dimethacrylates [5-20%] - Barium aluminosilicate glass [50-70%] - Amorphous silica [1-10%] - Camphorquinone [CQ] <1% - Accelerators <1% - Pigments <1% - BHT <0,1%

The total amount of inorganic filler is approximately 65 w% and the particle size ranges from 0.01 µm to 6 µm.

Clinical procedures

Trial fitting and checking the shade

1. Try the prosthodontic element in mouth to verify its fit on the prepared tooth.
2. Check the selected shade, with the Try In Paste. The shade of the Try In Paste matches the same shade of the cured INNOCEM Veneer LC composite.

Tooth structures preparation

1. Prepare the teeth in a conventional manner. Enamel surfaces should be etched, to improve the adhesive strength.
2. Try in restoration to verify its good adaptation, especially at the limits of the preparation.
3. Verify the choice of the shade of INNOCEM Veneer LC using the Try-In Paste.
4. Clean the tooth structures, rinse thoroughly with a water spray, and gently air dry. Do not over dry. The tooth structures should be just dry enough that the surface has a slightly glossy appearance.
5. Isolate the operating field (dental dam).
6. Apply a uniform coating of SEALBOND Universal onto the prepared dentin and enamel, using a disposable applicator for 15 seconds with a rubbing motion.
7. Gently dry with water-free and oil-free air for 5 to 10 seconds to evaporate the solvent base and to disperse the adhesive homogeneously. Light-cure 10 seconds.
8. For veneers and Maryland bridges, tooth structure must be etched.

Note
HEALBOND Universal provides high enamel bond strengths without the use of etching. However, etching gel can be used on enamel and/or dentin, especially in case of veneers with poor retention, and Maryland bridges. In that case, after step 1.3, etch Enamel or Dentin surfaces for 15 seconds, then rinse thoroughly ensuring that all etchant has been removed. Proceed to step 4.

Preparation of restorations

Treatment of Lithium Disilicate, Hybrid ceramic restorations

1. Etch the restoration with hydrofluoric acid according to the instructions of the manufacturer of the restorative material.
2. Rinse thoroughly with water for 15 seconds and dry with air, free of water and oil.
3. Apply a thin layer of ZC Primer to the adherent surface of the restoration using a micro-tip applicator for 15 seconds with a rubbing motion and allow the material to react for 60 seconds.

4. Dry, preferably with warm air (hand-dryer), until the solvent has evaporated. Do not rinse.

Treatment of Zirconia-based, Alumina-based and Composite restorations

1. Sandblast the internal surface of the restorations with 50µm alumina with a pressure of about 30 psi [0,2 MPa].
2. Rinse thoroughly with water for 15 seconds, if necessary, clean the restoration in an ultrasonic unit for about 1 minute and dry with air, free of water and oil.
3. Apply a thin layer of ZC Primer to the adherent surface of the restoration using a micro-tip applicator for 15 seconds with a rubbing motion and allow the material to react for 60 seconds.
4. Dry, preferably with warm air (hand-dryer), until the solvent has evaporated. Do not rinse.

Concomitant of prosthodontic element

1. Apply the chosen shade of INNOCEM Veneer LC on the prosthodontic element.
2. Gently place it on the tooth, leaving the paste overflow the whole margins.
3. Remove most of the excesses with a little disposable brush, holding the restorative element in place.
4. Residual excesses will still overflow. Light-cure through the prosthodontic element, using the "normal mode" of a LED curing light [1200 to 1400mW/cm²], for about 5 seconds.
5. After removing of the paste excesses (for example with a n°12 scalpel blade), light-cure the entire surface using "normal mode" of the LED curing light 1200-1400 mW/cm²), as well as the margins of the prosthodontic element, for 5 seconds/zone, on 3 to 4 zones,

depending of the size of the prosthodontic element (for example 4 x 5 seconds if the prosthodontic element is divided in 4 zones, with 4 second of pause between each high irradiance to light exposure, to let the tooth cooling down).

6. The use of a glycerine gel reduces the oxygen inhibition layer at the limits of the preparation.

Shaping

1. Shape any remaining uneven areas on the edges of prosthodontic elements and polish the marginal area with rubber cups filled with diamond polishing paste for ceramic. Use specific rubber cups for composites.
2. Check the occlusion.

Note

This device is used on young adults to elderly people, since definitive prosthodontic sessions are only conducted on adults.

Storage & shelf life

- The product must be used by the expiration date indicated on the package.
- Do not expose the product to light or warm sources.
- The product should be tightly closed immediately after use to prevent evaporation of volatile components.
- Store at a temperature between 5°C and 25°C. When the product is stored in the refrigerator, it should be brought back to room temperature 15 minutes minimum before use, in order to restore its normal viscosity and curing time.
- Use at room temperature
- Close immediately after use.

Important notice

- Do not throw away the box and the instruction for use until the product is finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

Dispose of the product in accordance with applicable local, regional and national regulations.

This material is intended exclusively for use in dentistry.

Processing must be done strictly in accordance with the instructions for use.

No liability can be accepted for any damage caused by the use of the instructions for use or the prescribed area of application. It is the responsibility of the user to check the material for its suitability and use for a purpose not expressly stated in the instructions for use.

Descriptions and information do not constitute a guarantee of properties and are not binding.

INNOCEM Veneer LC - Curing time				
Type and shade of the prosthodontic element	Maximum prosthodontic element thickness	Type of curing light and light intensity		
		"Boost" mode Blue LED & Blue & UV LED (≥ 2500 mW/cm²)	"Normal" mode Blue LED 1 Blue & UV LED (≥ 2500-1400 mW/cm²)	Halogen lamp (≥ 400 mW/cm²)
Veneers		Curing time		
Silica-based ceramic (lithium disilicate)	< 2 mm	5 sec/zone (3 or 4 zones, depending of the size of the prosthodontic element		
Hybrid ceramic, composite		10 sec		
Inlays, onlays		Curing time		
Silica-based ceramic, hybrid ceramic, composite resin	< 2 mm	5 sec/zone (3 or 4 zones, depending of the size of the prosthodontic element	15 sec	20 sec
Shades and Depth of cure		Bright, Light, Medium		

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