

# DUROC CLEAR

Translucent resin for surgical guide

Formulation without addition of TEGDMA, BisGMA, BisGMA



## ADVANTAGES

Resistant but slightly resilient, for an easy removal

A new layer can be applied on top of an existing one and stick to it

**Translucent** for a better vision of the surgical field

Adequate setting time

Low exothermicity



## INDICATIONS

Surgical guides directly in mouth

Stabilization of SICAT plates (Sirona)

Bite registration when unstable

## REFERENCES & PRESENTATIONS

**DUCL-50** 50ml cartridge + 10 mixing tips

**EJ-100** 100 yellow mixing tips

MD : Class IIa certified by Tüv Rheinland (0197), excl. tips, Class 1.

## MAKING OF A SURGICAL GUIDE DIRECTLY IN MOUTH IN A CAD/CAM ENVIRONMENT

